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Table (1): The types and sizes of sample containers, along with the procedures for sample storage and preservation.

Parameters/ matrix	Method/ Instrument	Container	Collected Volume	Preservation	Minimum sample size	Holding Time
Trace metals in water	ICP-MS	Teflon	500 ml	HNO ₃ , pH < 2 , fridge	50 ml	6 Months
Trace metals in sediment, soil	ICP-MS	Glass	2 Kg	Fridge	25 g	6 Months
Mercury in water	Mercury analyzer	Amber Glass	250 mL	HNO ₃ , pH < 2 , fridge	10 ml	28 Days

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Mercury in sediment, soil	Mercury analyzer	Amber Glass	250 mL	HNO ₃ , pH < 2 , fridge	10 ml	28 Days
Diesel Range Organics (DRO) in water	GC-FID	Amber Glass	1 L	HCL or H ₂ SO ₄ , pH < 2 , fridge	1L	14 Days
Polychlorinated biphenyls (PCBs) in water	GC-FID GC-MS/MS	Amber Glass	1 L	HCL, pH < 2 , fridge	1L	14 Days
Total Organic Carbon (TOC) in water	TOC Analyzer	Amber Glass	250 mL	HCL or H ₂ SO ₄ , pH < 2 , fridge	10 ml	28 Days
Polyaromatic hydrocarbon (PAH) in water	GC-MS	Amber Glass	1 L	HCL or H ₂ SO ₄ , pH < 2 , fridge	500 ml	14 Days
Oil & Grease in Water	Solid phase Extraction	Amber Glass (Brown)	2.5 L	HCL PH < 2 , fridge	1.5 L	28 Days
Total suspended solids (TSS) in Water	Filtration Classical Method	Glass	1 L	Fridge	1 L	Must be analyzed as soon as possible
pH	Probe	Glass	1 L	Fridge	250 mL	Must be analyzed as soon as possible (15 min.)
EC, TDS, & Salinity(Chloride) in Water	Probe	Glass	1 L	Fridge	250 mL	Must be analyzed as soon as possible (15 min.)
Residual Chlorine	Titration	Glass	1 L	Bottle should be filled (without air gap) tightly closed	250 mL	Must be analyzed as soon as possible (Analyze Immediately)
Chemical Oxygen Demand (COD) in Water	HACH Spect.	Amber Glass (Brown)	250 mL	H ₂ SO ₄ (0.5 mL/250 mL) PH < 2 , fridge	250 mL	28 Days

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Biochemical Oxygen Demand (BOD) in Water	BOD Probe	Amber Glass (Brown)	2.5 L	Incubator for 5 days at 20°C	1 L	Must be analyzed as soon as possible
Dissolved Oxygen (DO) in Water	Titration Classical Method	Glass	300 mL	2mL W1 + 2mL W2 at the field, , fridge	300 mL	Analyze Immediately
Total Kjeldahl Nitrogen (TKN) in Water	TKN Analyzer	Glass	1 L	H ₂ SO ₄ (0.5 mL/250 mL) PH < 2 , fridge	250 mL	28 Days
Ammonia (NH ₄ -N) in water	Spectrophotometric method	Glass	1 L	Filtration, Freezer	250 mL	Must be analyzed as soon as possible
Ortho-phosphate (PO ₄ -P) in Water	Spectrophotometric method	Glass	1 L	Filtration, Freezer	250 mL	Must be analyzed as soon as possible
Nitrite (NO ₂ -N) in Water	Spectrophotometric method	Glass	1 L	Filtration, Freezer	250 mL	Must be analyzed as soon as possible
Nitrate (NO ₃ -N) in Water	Spectrophotometric method	Glass	1 L	Filtration, Freezer	250 mL	Must be analyzed as soon as possible
Silicate (SiO ₂ -Si) in Water	Spectrophotometric method	Plastic	500 mL	Filtration, Freezer	250 mL	Must be analyzed as soon as possible
Total organic matter (TOM%) in Sediment, Soil	Titration	Glass	2 Kg	Freezer	500 gm	28 days
Total Calcium Carbonate (TCaCO ₃) % in Sediment, Soil	Titration	Glass	2 Kg	Freezer	500 gm	-
Grain size distributed (GR %) in Sediment, Soil	Sieving	Glass	2 Kg	Freezer	1 Kg	-
Chlorophyll a in Water	Spectrophotometric method	Plastic Black	2.5 L	1% MgCO ₃ Solution	1.3 L	Within (2 Days)

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Total & Fecal Coliform in Water	IDEXX method	Pre-sterilized plastic bottles with sod. thiosulphate	100 mL	At 4-10°C in ice	100 mL	6 hours of collection
Phytoplankton in Water	Compound Microscope	Plastic	1L	4% Formalin Solution	1L	12 months
Zooplankton in Water	Compound Microscope	Plastic	1L	4% Formalin Solution	1L	12 months
Benthic Flora in Sediment	By Stereo Microscope & Magnifying glass	Glass	2 Kg	Freezer	2 Kg	Keep in freezer
Benthic Fauna in Sediment	By Stereo Microscope & Magnifying glass	Glass	2 Kg	Freezer	2 Kg	Keep in freezer

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Table (2): Reference test methods used in the Environmental Laboratory Section.

No.	Test Code	Test Description	Test Method
Trace analysis test procedures		التحليل الكيميائي الدقيقة	
1.	SOP-TOC-Water-01	Determination of TOC in water samples by TOC-L analyzer	APHA -5310 B
2.	SOP-METALS-SOIL-02	Simultaneous Determination of Heavy Metals in Soil and Sediments Samples by ICP-MS	USEPA- 200.8
3.	SOP-METALS-WATER -03	Simultaneous Determination of Heavy Metals in water Samples by ICP-MS	USEPA- 200.8
4.	SOP-Hg-WATER-04	Determination of Mercury in water samples by CV AAS	USEPA 245.2
5.	SOP-HG-SED-05	Determination of Mercury in Sediment samples by CV AAS	USEPA 245.5

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6.	SOP-PAHS-WATER-06	Determination of PAHs in water samples by GC-	USEPA-8270
7.	SOP-PAHS-SED-07	Determination of PAHs in Sediment samples by GC-	USEPA-8270C
8.	SOP-TPHs-water-08	Determination of TPHs in water samples by GC-FID.	USEPA-8015B
9.	SOP-TPHs-SED-09	Determination of TPHs in Sediment samples by GC-FID.	UNEP-20
10.	SOP-PCBs-WATER-10	Determination of Polychlorinated Biphenyls (PCBs) in Water by Gas Chromatography Mass Spectrometry (GC-MS/MS)	ASTM-D8543
11.	SOPPCBs-SED-11	Determination of PCBs-Congeners in Sediment samples by GC-ECD.	USEPA-8082
12.	SOP-G&G-WATER-12	Determination of Oil & Grease in water samples by gravimetric method	USEPA 1664
التحاليل الفيزيوكيميائية Physio-chemical test procedures			
13.	SOP-DO-WATER-13	Determination of Dissolved Oxygen (DO) in water	APHA 4500 – O G
14.	SOP-BOD5-WATER-14	Determination of BOD5 in water samples by membrane electrode.	APHA 5210-B ISO 5815-1
15.	SOP-COD-WATER-15	Determination of COD in water samples by Closed Reflux spectrometric method	APHA 5220 D USEPA 8000 D
16.	SOP-COD- SALT WATER-16	LCK1814 COD cuvette test for salt water, 7-70 mg/L O ₂	ISO 6060-1989
17.	SOP-AMMONIA-WATER-17	Determination of Ammonia-N in water	USEPA 350.1 ISO 11732:1997
18.	SOP-NO ₃ -N-WATER-18	Determination of NO ₃ -N N-VCl ₃ -HCl-SW-Method	USEPA 353.1
19.	SOP-NO ₃ -N-WATER-19	Determination of NO ₃ -N Cadmium Reduction Method	USEPA 353.1 SM 4500 – NO ₃ G
20.	SOP-NO ₂ -N-WATER-20	Determination of NO ₂ -N in water	USEPA 354.1

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21.	SOP-PHOS-WATER-21	Determination of Phosphates in water	USEPA 365.1
22.	SOP-SILICATE-WATER-22	Determination of Silicates in water	USEPA 370.1
23.	SOP-TSS-WATER-23	Determination of Total Suspended Solids (TSS)	APHA-2540-D ISO 11923
24.	SOP-TDS-WATER-24	Determination of Total Dissolved Solids (TDS)	APHA 2510 B
25.	SOP-pH-WATER-25	Determination of pH in water samples	APHA 4500-H ⁺ B
26.	SOP-COND-WATER-26	Determination of electric conductivity	APHA 2510 B
27.	SOP- CHLORIDE - 27	Determination of chloride	APHA 4500 Cl-B
28.	SOP SALINITY - 28	Determination of Salinity	APHA 2520 B
29.	SOP-GRIN-SED-29	Determination of Grain Size in soil & sediment	ASTM - D422 - Standard Test Method for Particle-Size Analysis of Soils
30.	SOP- % (TCaCO ₃) -SED-30	Determination of % TCaCO ₃ in soil & sediment	Diagnosis and Improvement of Saline and Alkali Soil - Agriculture Handbook No. 60 -Richards 1954, (23b) p105. Titration method.
31.	SOP- TOM %-SED-31	Determination of % TOM in soil & sediment	Diagnosis and Improvement of Saline and Alkali Soil - Agriculture Handbook No. 60 -Richards 1954, (24) organic matter p105. Walkley – Black Wet Combustion Method
32.	SOP- Tar Balls -SED-32	Determination of Tar Balls (gm) in beach sediment	MOOPAM **, 3rd Edition, 1999
Biology and microbiology test procedures التحاليل البيولوجية			

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33.	SOP-T Bact-WATER-33	Determination of Total bacteria in water sample by IDEXX	APHA9223-B
34.	SOP-E COLI-WATER-34	Determination of Fecal Coliform as E.Coli in water sample by IDEXX	APHA9223-B
35.	SOP-PHYTO-WATER-35	Identification, Taxonomy & Numerical abundance of Phytoplankton Species	Phytoplankton Manual, Unesco, 1978.
36.	SOP-CHLORO-WATER-37	Spectrophotometric Determination of Chlorophyll-a in Water	ASTM-D373
37.	SOP-BENTIC-SED-38	Identification, Taxonomy & Numerical abundance of Benthic Fauna & Flora Species	Standard methods for the Examination of Water & Wastewater (SMWW) * 20th edition
38.	SOP- Enterococci -WATER-41	Detection of Enterococci by IDEXX Enterolert Quanti-Tray method	APHA 9230D
39.	SOP-MW-42	Microwave-based total acid digestion of soil and sediments	US EPA 3051A

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Table (3): International System of Units (SI) applied in reporting the results for seawater, wastewater, and groundwater analyses.

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Table (4): International System of Units (SI) applied in reporting the results of sediments.

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